6A53, 6A54

DE Gebrauchsanweisung	3
EN Instructions for Use	8
FR Instructions d'utilisation	14
IT Istruzioni per l'uso	20
ES Instrucciones de uso	26
PT Manual de utilização	33
NL Gebruiksaanwijzing	39
SV Bruksanvisning	45
DA Brugsanvisning	51
NO Bruksanvisning	57
FI Käyttöohje	62
PL Instrukcja użytkowania	68
HU Használati utasítás	74
CS Návod k použití	80
RO Instrucțiuni de utilizare	86
HR Upute za uporabu	92
SL Navodila za uporabo	97
SK Návod na používanie	103
BG Инструкция за употреба	109
TR Kullanma talimatı	115
EL Οδηγίες χρήσης	120
RU Руководство по применению	127
JA 取扱説明書	133
ZH 使用说明书	139
KO 사용 설명서	144

1 Produktbeschreibung

INFORMATION

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INFORMATION

Date of the last update: 2012-12-11

► Please read this document carefully.

► Follow the safety instructions.

1.1 Design and Function

The 6A53 and 6A54 Sliding Adapters permit the distal components of the modular prosthesis to be shifted in parallel, regardless of the angle adjustment by the pyramid adapters. The adjustments can be made either in the frontal plane (medial or lateral) or in the sagittal plane (anterior or posterior). The sliding adapters are equipped with different connections:

Product reference number	Connections
6A53	pyramid adapter/pyramid adapter
6A54	pyramid adapter/pyramid receiver

1.2 Combination Possibilities

For combination possibilities, please see the catalogue 646K2* or consult the manufacturer.

2 Application

2.1 Intended Use

The product is intended solely for lower limb prosthetic fittings.

2.2 Area of Application

Area of application according to the MOBIS® mobility system:

Approved for a body weight of up to **125 kg**.

The products may only be used in TT prostheses.

2.3 Environmental Conditions

Allowable Environmental Conditions
Temperature range for use: -10 °C to +60 °C (14 °F to 140 °F)
Allowable relative humidity 0 % to 90 %, non-condensing
Unallowable environmental conditions
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e. g. talcum)

2.4 Service Life

This component has been manufacturer-tested with 3 million load cycles according to ISO 10328. Depending on the level of user activity, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

⚠ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

3.2 General Safety Instructions

⚠ CAUTION

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with **"Inform the patient"**.

⚠ CAUTION

Excessive strain on the product

Fall due to breakage of load-bearing components

- ▶ Use the prosthesis components according to the MOBIS classification (see section "Area of Application").
- ▶ **Inform the patient.**

⚠ CAUTION

Unallowable combination of prosthesis components

Fall due to breakage or deformation of the product

- ▶ Only combine the product with prosthesis components that are approved according to the section "Combination Possibilities".
- ▶ Based on the instructions for use of the prosthesis components, verify that they may be combined with each other.

⚠ CAUTION**Exceeding the service life**

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded (see section "Service Life").
- ▶ **Inform the patient.**

NOTICE**Use under unallowable environmental conditions**

Damage to product due to unallowable environmental conditions.

- ▶ Do not expose the product to unallowable environmental conditions (see section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).
- ▶ **Inform the patient.**

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthesis components relative to each other and by the development of noises.

4 Scope of Delivery

The scope of delivery for the product is shown on page 2 (Figures 1 and 2).

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

6A53 Sliding Adapter					
Fig.		No.	Quantity	Designation	Reference number
–	■	–	1	Instructions for Use	647G763
1	–	–	1	Sliding Adapter	–

6A54 Sliding Adapter					
Fig.		No.	Quantity	Designation	Reference number
–		–	1	Instructions for Use	647G763
2	–	–	1	Sliding Adapter	–
2		–	4	Grub Screw	506G3=M8x12

5 Preparation for Use

CAUTION

Incorrect alignment or assembly

Injuries due to fall resulting from damaged prosthesis components

- Observe the alignment and assembly instructions.

5.1 Mounting the Sliding Adapter

CAUTION

Improper assembly of the screw connections

Fall due to breakage or loosening of the screw connections

- Clean the threads before every installation.
- Apply the installation torques provided.
- Observe the instructions for securing the screw connections and the use of the correct length.

Adjusting the Sliding Adapter and Mounting the Pyramid Adapter

A pyramid adapter adjustable in 1 mm increments is located on the side of the sliding adapter with the serrated surface. The pyramid adapter can be aligned proximally or distally as desired.

> Required tools and materials:

710D4 Torque Spanner, 636K13 Loctite® 241

- 1) Align the sliding adapter in the desired displacement direction (**a – p** or **m – l**).
- 2) Lift the ring on the adjustable pyramid adapter slightly until it no longer engages in the serrated surface (see Fig. 3).
- 3) Set the desired displacement and press the ring back down on to the serrated surface.
- 4) Connect the pyramid adapter with the applicable connecting component (see Fig. 4).

Mounting the Second Connection

6A53

> **Required tools and materials:**

710D4 Torque Spanner, 636K13 Loctite® 241

- Connect the pyramid adapter with the applicable connecting component.

6A54

> **Required tools and materials:**

710D4 Torque Spanner, 636K13 Loctite® 241

1) **For definitive mounting:**

Use Loctite® to secure the grub screws.

2) Screw in the grub screws.

3) **For trial fitting:**

Use the torque spanner to tighten the grub screws (tightening torque: **10 Nm**).

For definitive mounting:

Pre-tighten the grub screws with the torque spanner (tightening torque: **10 Nm**) and then tighten (tightening torque: **15 Nm**).

- 4) Replace any grub screws that are protruding or recessed too much with suitable ones (see selection table).

Selection Table for Set Screws	
Reference Number	Length (mm)
506G3=M8x10	10
506G3=M8x12-V	12
506G3=M8x14	14
506G3=M8x16	16

6 Maintenance

- The prosthesis components should be inspected after the first 30 days of use.
- Inspect the entire prosthesis for wear during normal consultations.
- Conduct annual safety inspections.

7 Legal Information

7.1 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with

tested modular components in accordance with the MOBIS Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorised by the manufacturer.

7.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

7.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction. Because warranty conditions may differ by country, please contact your local manufacturer distribution company for more information.

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Otto Bock HealthCare GmbH